

Delivery System

A-dec 300 Upgrade



A-dec 300 support center
and A-dec 300 Traditional
Delivery System.



Continental Delivery System



Standard Touchpad



Deluxe Touchpad

Your solution to conveniently deliver what you need, where you need it.

Small on the outside, big on the inside. With the A-dec 300 delivery system, you have a simple option to update your dental chair with genuine integration and easy positioning.

Add technology and improve ergonomics. Choose between Continental or Traditional. Either way, it's a delivery system designed to retrofit most dental chairs from leading manufacturers.

Your Authorized
A-dec Dealer is
ready to explain
how to upgrade
your existing chair
with the A-dec 300
Delivery System.



Seamless Integration

- Integrate up to two accessories, including intraoral cameras, electric handpieces, curing lights or ultrasonic instruments.
- Add future products easily with the A-dec 300's advanced delivery system platform.
- Control the cuspidor, dental light, handpieces, and integrated ancillaries with the optional deluxe touchpad.

Ergonomics & Reliability

- Reduce strain and fatigue. The A-dec 300 delivery system positions easily and places your handpieces and ancillaries within gentle reach.
- Get the reliability you've come to expect from a worldwide leader in quality dental solutions.

Options & Versatility

Also available on the A-dec 300 support center:

- **Assistant's holders**—choose 3- or 4-position holders with either cuspidor-mounted or telescoping arm.
- **Cuspidor**—the contemporary design features a removable contemporary-style, glass bowl.
- **Dental Light**—choose from A-dec's selection of post-mounted lights that meet your clinical needs, including the A-dec LED.
- **Monitor Mount**—positions easily over the patient for convenient patient consultation and education.

Learn More

For details, contact your local A-dec dealer. To request a list of authorized dealers, go to a-dec.com or call 1.800.547.1883.



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132
United States

Tel: 1.800.547.1883 within USA/CAN
Tel: 1.503.538.7478 outside USA/CAN
Fax: 1.503.538.0276
a-dec.com

A-dec Australia

Tel: 1.800.225.010 within AUS
Tel: +61.(0)2.8332.4000 outside AUS

A-dec China

Tel: +1.503.538.7478

A-dec United Kingdom

Tel: 0800.ADEC.UK (2332.85) within UK
Tel: +44.(0).24.7635.0901 outside UK

A-dec, the A-dec logo, Cascade, Century Plus, Continental, Decade, ICX, ICV, Performer, Preference, Preference Collection, Preference ICC, and Radius are trademarks of A-dec Inc. and are registered in the United States and other countries. A-dec 500, A-dec 400, A-dec 300, A-dec 200, A-dec Inspire, and reliablecreativesolutions are also trademarks of A-dec Inc. None of the trademarks or trade names in this document may be reproduced, copied, or manipulated in any manner without the express, written approval of the trademark owner.

Certain touchpad symbols are proprietary to A-dec Inc. Any use of these symbols, in whole or in part, without the express written consent of A-dec Inc., is strictly prohibited.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2017 A-dec Inc. All rights reserved. 85.6158.00/OP/1M/1-17/Rev D

